

Multi-Center Flexible Endoscope Channel Inspection Survey


59%

of Inspected Flexible Scope Working Channels Were Dirty

Participating Hospitals • 41 Endoscopes Inspected with a 2mm Waterproof Lumen Inspection System • 126

"Tools such as video boroscopes of an appropriate dimension (length and diameter) may be used to visually inspect the internal channels of some medical devices."

ANSI/AAMI ST91:2015 Flexible and semi-rigid endoscope processing in health care facilities, 12.4.2, Cleaning verification



"Visual inspection is recommended to make sure the endoscope is visibly clean."

Standards of Infection Prevention in Reprocessing of Flexible Gastrointestinal Endoscopes, 2016, Society of Gastroenterology Nurses and Associates (SGNA)



Number of Scopes Inspected (Percentage Dirty)

Duodenoscopes	48% (24 of 50)
Gastrosopes	65% (24 of 37)
Colonoscopes	67% (26 of 39)

Survey performed with INSPEKTOR® Channel Inspection System

SP Concepts www.inspektor.com

