

NEW 10" INSPEKTOR: 3X Larger Image



INSPEKTOR® CT 1.9 and DUO Video Borescope Inspection System




The Waterproof Video Borescope for use in GI and SPD Lumen Inspection

- A Plug-N-Play System with 10" monitor - Computer Not Required
- All-in-one Borescope Inspection System
- Waterproof Monitor and Borescope – For Use on Dirty or Clean Side
- Automatic White Balance and Illumination - Saves Time
- Small Enough to get into Bronchoscopes and Long Enough for GI Endoscopes



**New Mobile
Inspection Cart**



Dirty Arthroscopic Shaver

Debris Left in Endoscope

Dirty Duodenoscope

Borescope Sizes

INSPEKTOR CT ●
1.9 mm

INSPEKTOR CT MINI ●
1.2 mm



New 10-inch INSPEKTOR® CT Color HD Monitor with built in light-source

INSPEKTOR® System - Configurations

CT 1.9 System

INS7000-10	HD Color Monitor (10.1-inch diagonal) with built-in light source
INS8007*	1.9mm x 1.8m CT Digital Borescope
INS6005	Wall Mount Articulating Arm

CT DUO® System

INS7004	HD Color Monitor (10.1-inch diagonal) with built-in light source
INS8007*	1.9mm x 1.8m CT Digital Borescope
INS8008	1.2mm x 1.0m CT Digital Borescope
INS6005	Wall Mount Articulating Arm

INSPEKTOR CTC DUO® System - Non-Monitor Version

For use with a separate HD monitor or with a computer.

INS7005	Non-Monitor Controller with built-in light source
INS8007*	1.9mm x 1.8m CT Digital Borescope
INS8008	1.2mm x 1.0m CT Digital Borescope
INS6005	Wall Mount Articulating Arm



*Optional Borescope Length for all listed systems: INS8009 1.9mm x 1.0 meter